






Medical Devices - Symbols

Symbole	Title of Symbol	Description of Symbol
	Do not use if package is damaged.	Indicates a medical device that should not be used if the package has been damage or opened.
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.
	Do not re-use	Indicates a medical device that is intended for one use or for use on a single patient during a single procedure
	Contains or presence of natural rubber latex	Indicates the presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device which may cause allergic reactions.
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.